

[All Products](#) / [Optomechanics](#) / [Tube System](#) / [C, S, and T-Mount Components](#)

TECHSPEC® Achromat Prototyping System



- Imaging Lens System for Prototyping
- Ideal for Identifying Desired Magnification and f/#
- Compatible with Helicoid Barrel for Use With C-Mount Cameras

Common Specifications

General

Function:

Kit

Technical Information

Products

Title	Stock Number	Price	Buy
Complete Achromat Prototyping Kit	#54-666	₹84,353 Qty 6+ ₹74,296	Contact Us
Cell with Two Retainers and Apertures for Achromat Proto Kit	#54-719	₹19,951 Qty 6+ ₹17,597	10 to 12 days
Cell with Two Retainers for Achromat Prototyping Kit	#58-255	₹11,149 Qty 6+ ₹9,822	10 to 12 days
Aperture Kit for Achromat Prototyping Kit	#54-720	₹11,149 Qty 6+ ₹10,372	10 to 12 days



Copyright 2023, Edmund Optics India Private Limited, #267, Greystone Building, Second Floor, 6th Cross Rd, Binnamangala, Stage 1, Indiranagar, Bengaluru, Karnataka, India 560038

Phone: 1-800-363-1992 :

www.edmundoptics.com