

[See all 11 Products in Family](#)

## 1.75" x 2.75" Glassine Bags (150 Pieces)



1.75" x 2.75" Glassine Bag

Stock #17-183 **20+ In Stock**

1  MRP ₹1,408

Price inclusive of all taxes

**ADD TO CART**

### Volume Pricing

Qty 1+	₹1,408 each
Need More?	<a href="#">Request Quote</a>

### Product Downloads

#### General

Quantity:  
~150

#### Physical & Mechanical Properties

Dimensions (inches):  
1.75 x 2.75

#### Regulatory Compliance

Certificate of Conformance:

[View](#)

---

**Country of Origin:**

United States

**Imported By:**

Edmund Optics India Private Limited  
267, Greystone Building, Second Floor,  
6th Cross Rd, Binnamangala,  
Stage 1, Indiranagar, Bengaluru,  
Karnataka, India 560038  
Phone: +91- 80-6845 0000

---

## Product Details

- Glassine Bags, Jewel Boxes, and Cloth Pouches for General Optics
- Non-Contact Impact Cases for Laser and Precision Optics
- Ideal for Storing or Transporting Optics

Optics Packaging is used to safely store and protect optics against environmental or incidental damage when not in use. Glassine bags, cloth pouches, and jewel boxes are available for storing uncoated or coated optics including lenses, mirrors, and filters. Non-contact impact cases designed to hold optics only by their bevels to prevent contact with the optical surfaces are available for more sensitive optics such as laser optics. Optics Packaging is ideal for prolonging the life of optics when used in conjunction with [optics handling supplies](#).

**Note:** Non-contact impact cases are designed to hold optics with a specific size and thickness. Smaller optics will not be securely held which may result in damage.

---

;