

[See all 24 Products in Family](#)

## 20 x 20mm Glass Checkerboard Target NIST Certified



Opal Checkerboard Calibration Target



Stock #29-425 [CONTACT US](#)

- 1 + MRP ₹1,04,076

Price inclusive of all taxes

**ADD TO CART**

Volume Pricing	
Qty 1-4	₹1,04,076 each
Qty 5+	₹98,872 each
Need More?	<a href="#">Request Quote</a>

### Product Downloads

### General

**NIST Certification:**

Yes

### Physical & Mechanical Properties

**Pattern Size (mm):**

10 x 10 ±0.010

20 x 20 ± 0.10      **Dimensions (mm):**

---

1.52 ± 0.10      **Thickness (mm):**

---

0.2 x 0.2 ± 0.002      **Check Size (mm):**

## Optical Properties

**Coating:**  
Low Reflective Chrome (reflectivity <25% at 550nm)

**Substrate:** □  
Float Glass

## Regulatory Compliance

**Certificate of Conformance:**  
[View](#)

**Country of Origin:**  
United States

**Imported By:**  
Edmund Optics India Private Limited  
267, Greystone Building, Second Floor,  
6th Cross Rd, Binnamangala,  
Stage 1, Indiranagar, Bengaluru,  
Karnataka, India 560038  
Phone: +91- 80-6845 0000

## Product Details

- Chrome Pattern on Dimensionally Stable Substrates (Opal and Glass)
- Optimized for Reflective Illumination Calibration of Imaging Systems
- Various Sizes for Small or Large Field of Views

Checkerboard Calibration Targets are used in the calibration of 2D industrial and machine vision cameras to detect image distortion or sensor issues to ensure measurement accuracy. These Checkerboard targets feature white ivory and opal glass substrates that remain rigid, therefore not affecting the pattern by flexure. These Checkerboard targets minimize the effect of shadowing from back reflections, making them ideal for reflective illumination calibration of imaging systems with wide field of views. These targets are available in a range of sizes to complement cameras with small or large field of views.

**Note:** Targets with a specification of "NIST Certification: Yes" include a calibration certificate.

## Technical Information

