

**TECHSPEC® 20mm Diameter, C-Mount Thick Optic Mount**



Stock #54-626 **3 In Stock**

- 1 + MRP ₹6,003

● Price inclusive of all taxes

**ADD TO CART**

Volume Pricing	
Qty 1-4	₹6,003 each
Qty 5-9	₹5,337 each
Qty 10-25	₹4,954 each
Need More?	<a href="#">Request Quote</a>

Product Downloads

**General**

Fixed **Type:**

Optic Mount (Fixed) **Function:**

**Type of Optics:**

