

[See all 13 Products in Family](#)

20X Objective, CFI Plan Achromatic

See More by [Nikon](#)



Stock #75-355 NEW **1 In Stock**

MRP ₹37,489

● Price inclusive of all taxes

ADD TO CART

Volume Pricing

Qty 1+	₹37,489 each
Need More?	Request Quote

Product Downloads

General

Model Number:
MRL00202

Compatible Tube Lens Focal Length (mm):
Focal Length: 200mm

Type:
Microscope Objective

Style:
InfinityCorrected

Manufacturer:

Physical & Mechanical Properties1.1 **Field of View (mm):**58.67 **Length excluding Threads (mm):**27 **Maximum Diameter (mm):**85 **Weight (g):****Optical Properties**0.17 **Compatible Cover Glass Thickness (mm):**0.32 **Horizontal Field of View, 1/2" Sensor:**0.44 **Horizontal Field of View, 2/3" Sensor:**20X **Magnification:**0.40 **Numerical Aperture NA:**1.2 **Working Distance (mm):**22 **Field Number (mm):**59.87 **Parfocal Length (mm):**N/A **Immersion Liquid:****Sensor**2/3" **Maximum Sensor Format:****Threading & Mounting**M25 x 0.75 **Mounting Threads:****Regulatory Compliance**[View](#) **Certificate of Conformance:**Japan **Country of Origin:**

Imported By:
 Edmund Optics India Private Limited
 267, Greystone Building, Second Floor,
 6th Cross Rd, Binnamangala,
 Stage 1, Indiranagar, Bengaluru,
 Karnataka, India 560038
 Phone: +91- 80-6845 0000

Product Details

- Exceptional Flat-Field Imaging
- High Numerical Apertures and Oil Immersion Options Available
- Wide Magnification Range (1X to 100X)

Nikon's CFI Plan Achromat Objectives deliver exceptional flat-field imaging ensuring sharp, distortion-free clarity across the entire field of view, making these objectives ideal for both visual inspection and high-precision digital imaging. With high numerical apertures and specialized oil immersion options available, these objectives offer enhanced resolution and light-gathering capability for demanding high-magnification applications. Nikon CFI Plan Achromat Objectives are available in 1X up to 100X magnification, providing solutions for low-magnification overviews or detailed high-resolution imaging. Color corrected for the entire visible spectrum; these objectives are suitable for brightfield and fluorescence observation in routine lab work and photomicrography.

Technical Information

