

[See all 7 Products in Family](#)

2.5lb Test Instrument



2.5lb Test Instrument

Stock **#18-486** **1 In Stock**

⊖ 1 ⊕ ₹36,000

ADD TO CART

Volume Pricing

Qty 1+	₹36,000 each
Need More?	Request Quote

Product Downloads

General

Note:

Used for testing severe abrasion resistance.
Calibration of spring is NIST traceable.

Regulatory Compliance

Reach 224:

[Compliant](#)

Certificate of Conformance:

[View](#)

Country of Origin:

Product Details

- Test Optical Coatings for Abrasion and Adhesion Resistance per MIL Specifications
- Includes Instructions and Certifications of Conformance
- Replacement Erasers, Cheesecloth, and Other Components Available

Coating Testing Kits for Military Specifications (MIL-SPEC) feature all the tools and supplies needed to test optical coatings for abrasion and adhesion resistance to MIL-C-675, MIL-M-13508, MIL-F-48616, and MIL-PRF-13830B standards. Two instruments are included for testing abrasion resistance: a brass one with a 2 ½ lb spring to test severe abrasion resistance and an aluminum one with a 1 ¼ lb spring to test moderate abrasion resistance. The eraser inserts included with these two instruments are certified to conform to MIL-E-12397 and the calibration of the springs is NIST traceable. Cheesecloth meeting MIL-CCC-C-440 for moderate abrasion testing and tape meeting MIL-AA113 Rev. D for adhesion testing are also included. Coating Testing Kits for Military Specifications are an essential optical testing tool for those working in optics manufacturing, coating, or QA to ensure coated optics meet expected performance specifications.

;