

3" Length Bag, Cotton-Tipped Swabs



3" Length Bag, Cotton-Tipped Swabs (1 Bag)

Stock **#56-926** **6 In Stock**

-

1

+

₹518

ADD TO CART

Volume Pricing	
Qty 1+	₹518 each
Need More?	Request Quote

Product Downloads

General	
Type:	1 Bag
Number of Swabs:	100
Physical & Mechanical Properties	
Length (inches):	3.00
Regulatory Compliance	

Certificate of Conformance:

View

Country of Origin:

United States

Imported By:

Edmund Optics India Private Limited

Product Details

- Lint-Free
- 3" or 6" Wooden Stems
- Bags, Boxes, and Cases Available

The Cotton-Tipped Swab Applicators are ideal for optical cleanup applications. They are offered with either 3" or 6" wooden stems. These non-sterile swabs are great for providing versatility in a multitude of applications. Cotton-Tipped Swab Applicators are lint-free to be better functional in cleaning applications. Each bag contains 100 swabs, each box contains 1,000 swabs, and each case contains 10,000 swabs.

;