

[See all 5 Products in Family](#)

5X Objective, Nikon MTI LE Plan BF

See More by [Nikon](#)



5X Objective, Nikon MTI LE Plan BF

Stock #29-220 **2 In Stock**

1 MRP ₹73,278

1 Price inclusive of all taxes

ADD TO CART

Volume Pricing

Qty 1+	₹73,278 each
Need More?	Request Quote

Product Downloads

General

Model Number:
MJD00055

Compatible Tube Lens Focal Length (mm):
Focal Length: 200mm

Type:
Microscope Objective

Style:
Infinity Corrected

Nikon **Manufacturer:**

Optical Properties

5X **Magnification:**

0.10 **Numerical Aperture NA:**

27.16 **Depth of Field (µm):**

31.0 **Working Distance (mm):**

22 **Field Number (mm):**

60 **Parfocal Length (mm):**

Threading & Mounting

M25x0.75 **Mounting Threads:**

Regulatory Compliance

[View](#) **Certificate of Conformance:**

China **Country of Origin:**

Imported By:
Edmund Optics India Private Limited
267, Greystone Building, Second Floor,
6th Cross Rd, Binnamangala,
Stage 1, Indiranagar, Bengaluru,
Karnataka, India 560038
Phone: +91- 80-6845 0000

Product Details

- High NA
- Long Working Distance
- Parfocal Distance of 60mm For Clear Images

Nikon CF160 LE Plan EPI Brightfield Infinity Corrected Objectives are designed for Nikon's CF160 system, a proprietary optical system combining infinity optics with the superior performance of the CF optical system. Featuring a parfocal distance of 60mm, these objectives have long working distances while maintaining a high numerical aperture, and producing clear images with high contrast and resolution. Various magnification options are available, ranging from 5X to 50X. Nikon CF160 LE Plan EPI Brightfield Infinity Corrected Objectives are ideal for industrial microscopy applications such as technical cleanliness, steel quality rating, and metallography.