

Arri Type PL Mount Adapter



Arri Type PL Mount Adapter

Stock **#16-141** **1 In Stock**

-

1

+

₹17,550

ADD TO CART

Volume Pricing	
Qty 1+	₹17,550 each
Need More?	Request Quote

Product Downloads

General

Note:

Converts T-Mount thread to ARRI PL Bayonet Mount with locking screws.

Regulatory Compliance

Compliant

RoHS 2015:

Compliant

Reach 224:

Certificate of Conformance:

View	
United States	Country of Origin:
Edmund Optics India Private Limited	Imported By:

- Designed for Use from Macro to Infinity
- Ideal for Industrial and Cinema Production Applications
- Excellent Sharpness Over Depth of Field
- **MikroMak™ PetzVex™ Prime Probe Lenses** Also Available

MikroMak Prime Probe Lenses are available in five different focal lengths: 40, 50, 70, 90, and 125mm. The 40 and 50mm lenses provide wide field images with moderate macro magnifications up to 1.75X (50mm) or 2X (40mm). The 40mm lens can be converted to have a 25mm focal length with 1X magnification through an adapter (sold separately). The 70, 90, and 125mm lenses provide high magnification and resolutions typical of objectives used on standard laboratory microscopes with a maximum magnification of 5.5X (70 and 125mm) or 6X (90mm). A suite of accessories including adapters to C-Mount, L-Mount, and Sony mounts, mounting solutions, spacers, and gears for motorization are available; see each specific accessory for more information.