

TECHSPEC®

Cell with Two Retainers and Apertures for Achromat Proto Kit



Stock #54-719

CONTACT US

-

1

+

₹23,542

Price inclusive of all taxes

ADD TO CART

Volume Pricing	
Qty 1-5	₹23,542 each
Qty 6-10	₹20,764 each
Qty 11-25	₹19,146 each
Need More?	Request Quote

Product Downloads

General	
Prototyping Kit	Type:
Kit	Function:

Regulatory Compliance	
RoHS 2015:	

