


[See all 21 Products in Family](#)

## Clear, 3cc Barrel (100 count), Dispenser Barrel



Stock #25-095 **2 In Stock**

⊖ 1 ⊕ MRP ₹17,735

● Price inclusive of all taxes

**ADD TO CART**

### Volume Pricing

Qty 1+	₹17,735 each
Need More?	<a href="#">Request Quote</a>

### Product Downloads

### Regulatory Compliance

Certificate of Conformance:

[View](#)

Country of Origin:

United States

Imported By:

Edmund Optics India Private Limited  
267, Greystone Building, Second Floor,  
6th Cross Rd, Binnamangala,

## Product Details

- Empty Barrels are Available in Clear, Amber, and Opaque Versions
- Tips are Available in Plastic and Steel
- [Dispensing Guns](#) and [Pre-Loaded Norland Adhesive Barrels](#) are Also Available

Dispenser barrels, tips, and other accessories are available for use with dispensing guns and hand plungers. Dispenser Barrels are available in a clear version for normal fluids, an amber version for UV-curing fluids, and an opaque version for light sensitive fluids. Dispenser tips are available in plastic and steel to ensure varying levels of controllable deposits without leakage.

**Note:** Barrels, tips, and adhesive are sold separately from [dispensers](#). [Pre-loaded Norland adhesive barrels](#) are also available.

---

;