

[See all 48 Products in Family](#)

TECHSPEC® Filter Cube for Olympus BX3/IX3 Microscopes



Olympus Filter Cube #86-833

Stock **#86-833** **20+ In Stock**

⊖ 1 ⊕ MRP ₹75,933

📌 Price inclusive of all taxes

ADD TO CART

Volume Pricing	
Qty 1+	₹75,933 each
Need More?	Request Quote

Product Downloads

General

Model Number:
N2711800

Type:
Fluorescence Filter Kit

Manufacturer:
Olympus

Dimensions of Compatible Optics (mm):
Excitation Filter: 25mm Dia.
Emission Filter: 25mm Dia.
Dichroic Mirror: up to 38 x26mm

Note:
Compatible with Olympus BX3 and IX3 Series
Microscope

Physical & Mechanical Properties

Dimensions (mm):
48 x 41 x 39

Threading & Mounting

Thickness of Compatible Optics (mm):
Excitation Filters: up to 6mm
Emission Filters: up to 6mm
Dichroic Mirror: up to 1.4mm

Regulatory Compliance

Certificate of Conformance:
[View](#)

Country of Origin:
Japan

Imported By:
Edmund Optics India Private Limited
267, Greystone Building, Second Floor,
6th Cross Rd, Binnamangala,
Stage 1, Indiranagar, Bengaluru,
Karnataka, India 560038
Phone: +91- 80-6845 0000

Product Details

- Ensure Proper Alignment of Filter Sets
- Easily Adapt to Common Microscopes
- Available in a Variety of Wavelengths
- Unmounted [TECHSPEC® Fluorescence Filter Kits](#) Also Available

TECHSPEC® Pre-Mounted Fluorescence Filter Cube Sets are designed to easily integrate into popular Olympus, Nikon, Zeiss, and Zaber microscopes. Each fluorescence filter cube set consists of one [TECHSPEC Fluorescence Dichroic Filter](#) and two [TECHSPEC Fluorescence Bandpass Filters](#) that have been specifically oriented to achieve optimum excitation and emission of popular fluorophores. Additionally these filters are also used in polymerase chain reaction (PCR) diagnostic instruments. For increased adaptability, filters can be easily removed and replaced to meet the needs of the application. Unmounted [TECHSPEC Fluorescence Filter Kits](#) are also available.

Built-to-Order Product: Each product is assembled from stock components at the time of your order. If all components are in-stock, your product will typically ship within 3 - 5 days. If you have any questions about delivery or inventory status, please [contact us](#).

Built-to-Order Product: Each product is assembled from stock components at the time of your order. If all components are in-stock, your product will typically ship within 3 - 5 days. If you have any questions about delivery or inventory status, please [contact us](#).