

**TECHSPEC® 1/4-20 Stud, English Post Kit**



Stock **#68-230** 10 TO 12 DAYS

- 1 +

**₹6,047.00**

ADD TO CART

|           |
|-----------|
| Qty 1+    |
| ₹6,047.00 |

|                               |
|-------------------------------|
| Volume Pricing                |
| <a href="#">Request Quote</a> |

|                       |
|-----------------------|
| Get Product Downloads |
|-----------------------|

**Regulatory Compliance**

RoHS 2015:  
[Compliant](#)

Certificate of Conformance:  
[View](#)