

High Performance Vibration Isolation System

See More by [TMC](#)



Stock **#87-226** [CONTACT US](#)

- 1 +

 ₹3,57,750

ADD TO CART

Volume Pricing	
Qty 1-4	₹3,57,750 each
Qty 5+	₹3,39,750 each
Need More?	Request Quote

Product Downloads

General

Type:

Universal

Physical & Mechanical Properties

Thickness (inches):

3.0

Length (inches):

20.0	
18.0	Width (inches):
Construction: Electropolished stainless steel sides and base Textured powder coat carbon steel top plate	
80 (nitrogen or air)	Feed Pressure (psi):
3.0	Height (inches):
67 (150 lb)	Load Capacity (kg):
22 (50 lb)	Weight (kg):

Regulatory Compliance	
Compliant	RoHS 2015:
View	Certificate of Conformance:
United States	Country of Origin:
Edmund Optics India Private Limited	Imported By:

Product Details

- Passive Vibration Isolation
- Lightweight, Compact Profile for Easy Portability
- Self-Leveling Design with Height Control Valves

The TMC™ High Performance Vibration Isolation System is designed to remove the unwanted effects of vibration from challenging microscopy applications. The lightweight, low-profile design weighs only 22kg, easing portability and lab integration. The TMC™ High Performance Vibration Isolation System consists of an electropolished stainless steel chassis and powder finish steel top plate. The self-leveling TMC™ High Performance Vibration Isolation System also includes height control valves for modifying the height of the platform.