

[See all 48 Products in Family](#)

TECHSPEC[®] ICG Filter Cube Set - Olympus



Olympus Filter Cube #86-833



Stock #86-381-OLY

CONTACT US

-

1

+

₹1,07,550

ADD TO CART

Volume Pricing	
Qty 1-9	₹1,07,550 each
Qty 10-25	₹97,830 each
Qty 26-49	₹92,340 each
Need More?	Request Quote

Product Downloads

General	
Fluorescence Filter Kit	Type:
Indocyanine Green (ICG)	Compatible Fluorophore:
#84-105	Excitation Filter:

#84-107	Emission Filter:
#86-336	Dichroic Filter:
Olympus	Compatible Microscope Brand:

Optical Properties

Fused Silica (Corning 7980)	Substrate: □
Hard Coated	Coating:
748 - 789	Excitation Wavelength (nm):
814 - 851	Emission Wavelength (nm):
801.00	Dichroic Cut-On Wavelength (nm):
748 - 851	Wavelength Range (nm):

Regulatory Compliance

View	Certificate of Conformance:
United States	Country of Origin:
Edmund Optics India Private Limited	Imported By:

Product Details

- Ensure Proper Alignment of Filter Sets
 - Easily Adapt to Common Microscopes
 - Available in a Variety of Wavelengths
 - Unmounted [TECHSPEC® Fluorescence Filter Kits](#) Also Available
- TECHSPEC® Pre-Mounted Fluorescence Filter Cube Sets are designed to easily integrate into popular Olympus, Nikon, Zeiss, and Zaber microscopes. Each fluorescence filter cube set consists of one [TECHSPEC Fluorescence Dichroic Filter](#) and two [TECHSPEC Fluorescence Bandpass Filters](#) that have been specifically oriented to achieve optimum excitation and emission of popular fluorophores. Additionally these filters are also used in polymerase chain reaction (PCR) diagnostic instruments. For increased adaptability, filters can be easily removed and replaced to meet the needs of the application. Unmounted [TECHSPEC Fluorescence Filter Kits](#) are also available.
- Built-to-Order Product: Each product is assembled from stock components at the time of your order. If all components are in-stock, your product will typically ship within 3 - 5 days. If you have any questions about delivery or inventory status, please [contact us](#).
- Built-to-Order Product:** Each product is assembled from stock components at the time of your order. If all components are in-stock, your product will typically ship within 3 - 5 days. If you have any questions about delivery or inventory status, please [contact us](#).