

[See all 15 Products in Family](#)

InFocus Module for Nikon/Olympus/Mitutoyo/Achroid

See More by [Infinity Photo-Optical Company](#)



InFocus Module for Nikon/Olympus/Mitutoyo/Achroid

Stock **#33-137** CLEARANCE CONTACT US

⊖ 1 ⊕ MRP ₹2,32,047

i Price inclusive of all taxes

ADD TO CART

Volume Pricing	
Qty 1+	₹2,32,047 each
Need More?	Request Quote

Product Downloads

General

Fixed Magnification Lens **Type:**

Nikon/Olympus/Mitutoyo/Achroid **Compatible Objectives:**

Physical & Mechanical Properties

73.00 **Length (mm):**

Regulatory Compliance

Compliant **RoHS 2015:**

Compliant **Reach 224:**

View **Certificate of Conformance:**

Country of Origin:
United States

Imported By:
Edmund Optics India Private Limited
267, Greystone Building, Second Floor,
6th Cross Rd, Binnamangala,
Stage 1, Indiranagar, Bengaluru,
Karnataka, India 560038
Phone: +91-80-6845 0000

Product Details

- Enable Dynamic Optical Focusing and Reduce Spherical Aberration
- Compatible With All Infinity Corrected Objectives
- Color Matched for Major Brand Objectives

InFocus™ Dynamic Optical Focusing Systems are used with microscopes to enable dynamic optical focusing, and drastically reduce the spherical aberration caused by the sample's cover glass thickness or suspension fluid. These dynamic optical focusing systems eliminate the need for correction collars and expensive motorized computer activation. InFocus™ Dynamic Optical Focusing Systems allows micromanipulation of the microscope, and can be integrated with any infinity corrected objectives. Variations in magnification are minimal over the range of focal translation. InFocus™ Dynamic Optical Focusing System modules are color corrected for Mitutoyo, Achrovid, Zeiss, Leica, Nikon, or Olympus objectives and can be connected to most major brand microscopes using a series of stand adapters (sold separately).

Technical Information

#33-137 &
#33-138

