

TECHSPEC® Mounting Clamp for 1 - 3X Research-Grade Variable Beam Expander



Mounting Clamp for 1 - 3X Research-Grade Variable Beam Expander, #89-294

Stock **#89-294** 10 TO 12 DAYS

⊖ 1 ⊕ **₹12,128.00**

ADD TO CART

Qty 1-4
₹12,128.00
Qty 5+
₹10,723.00
Volume Pricing
Request Quote

Get Product Downloads



Type:
Mounting Clamp

Regulatory Compliance

RoHS:
[Compliant](#)

