

[See all 15 Products in Family](#)

# Nikon CF160 Objective Adapter (M25)

See More by [Infinity Photo-Optical Company](#)



InFocus Module for Leica, #33-139



Stock #54-589 [CONTACT US](#)

1 MRP ₹10,301

● Price inclusive of all taxes

**ADD TO CART**

## Volume Pricing

Qty 1+	₹10,301 each
Need More?	<a href="#">Request Quote</a>

## Product Downloads

## General

Type: Lens Accessory

## Regulatory Compliance

Compliant

RoHS 2015:

## Certificate of Conformance:

## Country of Origin:

United States

## Imported By:

Edmund Optics India Private Limited  
 267, Greystone Building, Second Floor,  
 6th Cross Rd, Binnamangala,  
 Stage 1, Indiranagar, Bengaluru,  
 Karnataka, India 560038  
 Phone: +91- 80-6845 0000

## Product Details

- Enable Dynamic Optical Focusing and Reduce Spherical Aberration
- Compatible With All Infinity Corrected Objectives
- Color Matched for Major Brand Objectives

InFocus™ Dynamic Optical Focusing Systems are used with microscopes to enable dynamic optical focusing, and drastically reduce the spherical aberration caused by the sample's cover glass thickness or suspension fluid. These dynamic optical focusing systems eliminate the need for correction collars and expensive motorized computer activation. InFocus™ Dynamic Optical Focusing Systems allows micromanipulation of the microscope, and can be integrated with any infinity corrected objectives. Variations in magnification are minimal over the range of focal translation. InFocus™ Dynamic Optical Focusing System modules are color corrected for Mitutoyo, Achrovid, Zeiss, Leica, Nikon, or Olympus objectives and can be connected to most major brand microscopes using a series of stand adapters (sold separately).

## Technical Information

