


[See all 4 Products in Family](#)

Olympus PLN 10X, Ph1, Phase Contrast Objective



Stock #86-935 **NEW** 1 In Stock

- 1 + ₹48,720

ADD TO CART

Volume Pricing	
Qty 1+	₹48,720 each
Need More?	Request Quote

Product Downloads

General

Model Number:
PLN10XPH

Type:
Microscope Objective

Style:
Phase Contrast

Manufacturer:
Olympus

Note:
Requires Ph1 phase annulus for phase contrast observation

Physical & Mechanical Properties

2.2 Field of View (mm):

34.27 Length excluding Threads (mm):

24 Maximum Diameter (mm):

Optical Properties

N/A Compatible Cover Glass Thickness (mm):

18.00 Focal Length FL (mm):

10X Magnification:

0.25 Numerical Aperture NA:

1.34 Resolving Power (μm):

4.40 Depth of Field (μm):

10.6 Working Distance (mm):

22 Field Number (mm):

45 Parfocal Length (mm):

N/A Immersion Liquid:

440 Depth of Focus (μm):

Threading & Mounting

RMS / 20.32mm x 36 TPI Mounting Threads:

Regulatory Compliance

[View](#) Certificate of Conformance:

Japan Country of Origin:

Edmund Optics India Private Limited Imported By:

Product Details

- Flat-Field Plan Achromat Design for Sharp, Distortion-Free Images
- Precision Phase Ring Matching (Ph1, Ph2, Ph3) For Optimal Contrast
- Available In Multiple Magnifications (10X, 20X, 40X) Plus A 100X Oil Immersion Option
- Optimized For the UIS2 Optical System, Ensuring Seamless Integration with Olympus CX43

Olympus PLN Phase Contrast Objectives provide high-performance imaging for precise cellular analysis on the [Olympus CX43](#) microscope. These objectives are engineered for accurate phase contrast imaging across magnifications, offering excellent contrast and flat-field correction. Designed specifically for the UIS2 optical system, they deliver sharp, high-clarity images ideal for live, unstained specimens. Olympus PLN Phase Contrast Objectives, such as the 10X Ph1 ([#86-935](#)) and 100X Ph3 oil ([#87-013](#)), ensure reliable results across a wide range of biological applications. Their consistent quality and compatibility with standard Olympus phase contrast condensers make them a trusted choice for professional labs.