

[See all 14 Products in Family](#)

# Olympus Super Wide Tube Lens

See More by [Olympus](#)



#23-723: Olympus Super Wide Tube Lens

Stock **#23-723** **1 In Stock**

1  MRP ₹99,828

**i** Price inclusive of all taxes

**ADD TO CART**

Volume Pricing	
Qty 1+	₹99,828 each
Need More?	<a href="#">Request Quote</a>

## Product Downloads

### General

N6388800 **Model Number:**

Tube Lens **Type:**

Olympus **Manufacturer:**

### Physical & Mechanical Properties

**Weight (g):**

## Optical Properties

180.00 **Focal Length FL (mm):**

1X **Magnification:**

400 - 700 **Wavelength Range (nm):**

## Regulatory Compliance

[View](#) **Certificate of Conformance:**

China **Country of Origin:**

**Imported By:**  
Edmund Optics India Private Limited  
267, Greystone Building, Second Floor,  
6th Cross Rd, Binnamangala,  
Stage 1, Indiranagar, Bengaluru,  
Karnataka, India 560038  
Phone: +91- 80-6845 0000

## Product Details

- Empty Filter Cube and [Fluorescence Filter Pre-Mounted Version](#) Available
- Wide Field and Single Port Tube Lenses
- Low Auto Fluorescence Immersion Oil for High NA [Objectives](#)
- C-Mount, T-Mount Adapters and Camera Coupler for System Integration

Olympus Tube Lenses and Accessories are specifically designed for use with Olympus microscopes or objectives.

The Olympus Empty Filter Cubes can accept 25mm diameter [excitation and emission](#) filters and 25.2 x 35.6mm [dichroic filters](#) for customized fluorescence configurations. The [#86-832](#) Filter Cube is compatible with Olympus BX2 (BX61/51/41) and IX2 (IX81/71/51) series microscopes. The [#86-833](#) Filter Cube is compatible with Olympus BX3 (BX63/53/43) and IX3 (IX83/73/53) series microscopes. [Pre-Mounted Fluorescence Filter Cube Sets](#) utilizing the [#86-833](#) Filter Cube are available for a variety of standard fluorophores.

Olympus Tube Lenses feature a 180mm focal length. The Olympus Wide Field Tube Lens ([#36-401](#)) offers a wide 26.5mm diameter image circle, and can be easily integrated into T-Mount systems with the Olympus Wide Field Tube Lens to T-Mount Holder ([#11-149](#))\*. The Olympus Single Port Tube Lens ([#86-835](#)) features a smaller 22mm diameter image circle, but is designed to easily integrate with microscope observation tubes for video observation. A C-Mount video microscope can be configured by using the [#86-835](#) tube lens, [#89-606](#) Olympus C-Mount Adapter, and [#89-605](#) Olympus 1X Camera Coupler.

Olympus Low Auto Fluorescence Immersion Oil [#86-834](#) is designed for use with Oil Immersion Objectives to achieve a high numerical aperture. The immersion oil is specially formulated for reducing autofluorescence that is typically associated with oils, making it an ideal choice for Fluorescence Microscopy applications.

\***Note:** The four socket head cap screws must be removed to mount the [#36-401](#) Olympus Wide Field Tube Lens in the [#11-149](#) Olympus Wide Field Tube Lens to T-Mount Holder. See Technical Information Tab.

## Technical Information

