


[See all 1 Products in Family](#)

Edmund Optics® Lab Storage Kit



Stock **#29-092** **5 In Stock**

1 MRP ₹9,558

Price inclusive of all taxes

ADD TO CART

Volume Pricing

Qty 1+	₹9,558 each
Need More?	Request Quote

Product Downloads

General

Contents of Kit:

24 component storage drawers
12 optics storage trays
drawer labels
marker

Regulatory Compliance

Certificate of Conformance:

[View](#)

United States **Country of Origin:**

Imported By:
Edmund Optics India Private Limited
267, Greystone Building, Second Floor,
6th Cross Rd, Binnamangala,
Stage 1, Indiranagar, Bengaluru,
Karnataka, India 560038
Phone: +91- 80-6845 0000

Product Details

- Maintain an Organized Optics Lab with Wall Mountable Storage
- 24 Component Storage Drawers, 12 Optics Storage Trays, Drawer Labels, and a Marker
- Convenient Storage for Lab Supplies
- Durable Boxes with Labels for Organization

Edmund Optics® Lab Storage Kit provide an ideal solution for protecting and organizing components in an optics lab environment. Each Lab Storage Kit includes 24 component storage drawers, 12 optics storage trays, drawer labels, and a marker. These convenient storage cases can be placed on a lab bench, easily wall mounted, and will help alleviate clutter and prevent loss and damage to precision components

;