

[See all 5 Products in Family](#)

Replacement Base for Peak Illuminated Direct Microscopes



Stock **#56-496** **1 In Stock**

- 1 + ₹3,848

ADD TO CART

Volume Pricing

Qty 1-4	₹3,848 each
Qty 5+	₹3,645 each
Need More?	Request Quote

Product Downloads

Regulatory Compliance

[View](#) Certificate of Conformance:

Country of Origin:
Japan

Imported By:
Edmund Optics India Private Limited

Product Details

- Multiple Magnification Options
- Illuminated Models Available with English or Metric Reticles
- Removable Clear Acrylic Base

Peak Direct Measuring Microscopes are high power, portable microscopes designed for precise depth measurements. Peak Direct Measuring Microscopes feature a 10X focusable eyepiece and a scale with an accuracy of 0.1mm for measuring the length, width, or depth of an object. The clear cutaway acrylic base allows for manipulation of the sample while viewing. Illuminated models feature a removable pen illuminator, and include two AA batteries. Peak Direct Measuring Microscopes are sold with a hard-shell vinyl storage case.

Note: Reticle should not be removed or replaced.
