

TECHSPEC[®] TRITC Filter Cube Set - Olympus



Olympus Filter Cube #86-833



Stock **#67-007-OLY** 10 TO 12 DAYS

-

1

+

₹1,07,550

ADD TO CART

Qty 1-9
₹1,07,550
Qty 10+
₹97,830
Volume Pricing
Request Quote

Product Downloads

SPECIFICATIONS

General

Fluorescence Filter Kit	Type:
TRITC	Compatible Fluorophore:
#67-032	Excitation Filter:
#67-034	Emission Filter:
#67-082	Dichroic Filter:
Olympus	Compatible Microscope Brand:
Optical Properties	
Fused Silica (Corning 7980)	Substrate: ☐
Hard Coated	Coating:
532 - 554	Excitation Wavelength (nm):
570 - 613	Emission Wavelength (nm):
562.00	Dichroic Cut-On Wavelength (nm):
532 - 613	Wavelength Range (nm):
Regulatory Compliance	
Not Compliant	RoHS:
View	Certificate of Conformance:

PRODUCT DETAILS

- Ensure Proper Alignment of Filter Sets
- Easily Adapt to Common Microscopes
- Available in a Variety of Wavelengths
- Unmounted [TECHSPEC® Fluorescence Filter Kits](#) Also Available

TECHSPEC® Pre-Mounted Fluorescence Filter Cube Sets are designed to easily integrate into popular Olympus or Nikon microscopes. Each fluorescence filter cube set consists of one [TECHSPEC Fluorescence Dichroic Filter](#) and two [TECHSPEC Fluorescence Bandpass Filters](#) that have been specifically oriented to achieve optimum excitation and emission of popular fluorophores. Additionally these filters are also used in polymerase chain reaction (PCR) diagnostic instruments. For increased adaptability, filters can be easily removed and replaced to meet the needs of the application. Unmounted [TECHSPEC Fluorescence Filter Kits](#) are also available.

Built-to-Order Product: Each product is assembled from stock components at the time of your order. If all components are in-stock, your product will typically ship within 3 - 5 days. If you have any questions about delivery or inventory status, please [contact us](#).