

Tube Lens



Stock #18-505 **1 In Stock**

⊖ 1 ⊕ ₹46,800

ADD TO CART

Volume Pricing	
Qty 1-4	₹46,800 each
Qty 5+	₹44,460 each
Need More?	Request Quote

Product Downloads

SPECIFICATIONS

General

Note:

One required between camera tube and tube lens spacer/filter cube

Physical & Mechanical Properties

Length (mm):
16.15

Maximum Diameter (mm):
36

Optical Properties

Focal Length FL (mm):
120.00

Magnification:
1X

Threading & Mounting

Mounting Threads:
M29 x 0.75 & Ø30mm Dovetail

Regulatory Compliance

Certificate of Conformance:
[View](#)

PRODUCT DETAILS


- Up to 43% Smaller than Industry Standard Systems
- Optimized Performance up to 1.1" Sensor
- Ideal for Life Science and Industrial Applications
- Winner of the [2023 Bronze VSD Innovators Award](#)

Edmund Optics® 120i Plan APO Infinity Corrected Objectives are infinite conjugate microscope objectives that offer exception image flatness and chromatic correction over 400-700nm. These objectives were designed for a 1.1" sensor to offer a large field of view and reliable image quality for high-throughput imaging applications. Ideal for life science and machine vision applications such as microscopy, flow cytometry, pharmaceutical inspection, and assembly line or fault inspection.

These objectives won the [2023 Bronze VSD Innovators Award](#).

TECHNICAL INFORMATION

Compact Tube Lens System



Standard Tube Lens System

