

[See all 22 Products in Family](#)

Tube Lens Spacer Assembly



Stock **#18-506** **1 In Stock**

1 ₹9,000

ADD TO CART

Volume Pricing	
Qty 1-4	₹9,000 each
Qty 5+	₹8,550 each
Need More?	Request Quote

Product Downloads

General

Note:

One required between tube lens and objective

Threading & Mounting

Mounting Threads:

C-Mount & M29 x0.75

Regulatory Compliance

Certificate of Conformance:

[View](#)

Country of Origin:

China

Imported By:

Edmund Optics India Private Limited

Product Details


- Up to 43% Smaller than Industry Standard Systems
- Optimized Performance up to 1.1" Sensor
- Ideal for Life Science and Industrial Applications
- Winner of the [2023 Bronze VSD Innovators Award](#)

Edmund Optics® 120i Plan APO Infinity Corrected Objectives are infinite conjugate microscope objectives that offer exceptional image flatness and chromatic correction over 400-700nm. These objectives were designed for a 1.1" sensor to offer a large field of view and reliable image quality for high-throughput imaging applications. Ideal for life science and machine vision applications such as microscopy, flow cytometry, pharmaceutical inspection, and assembly line or fault inspection.

These objectives won the [2023 Bronze VSD Innovators Award](#).

Technical Information

Compact Tube Lens System



Standard Tube Lens System

