

[See all 22 Products in Family](#)

Tube Lens



Stock #18-505 **3 In Stock**

MRP ₹55,224

Price inclusive of all taxes

ADD TO CART

Volume Pricing	
Qty 1-4	₹55,224 each
Qty 5+	₹52,463 each
Need More?	Request Quote

Product Downloads

General

Note:

One required between camera tube and tube lens spacer/filter cube

Physical & Mechanical Properties

Length (mm): 16.15

Optical Properties

120.00 Focal Length FL (mm):

1X Magnification:

Threading & Mounting

M29 x 0.75 & Ø30mm Dovetail Mounting Threads:

Regulatory Compliance

[View](#) Certificate of Conformance:

China Country of Origin:

Imported By:
Edmund Optics India Private Limited
267, Greystone Building, Second Floor,
6th Cross Rd, Binnamangala,
Stage 1, Indiranagar, Bengaluru,
Karnataka, India 560038
Phone: +91-80-6845 0000

Product Details


- Up to 43% Smaller than Industry Standard Systems
- Optimized Performance up to 1.1" Sensor
- Ideal for Life Science and Industrial Applications
- Winner of the [2023 Bronze VSD Innovators Award](#)

Edmund Optics® 120i Plan APO Infinity Corrected Objectives are infinite conjugate microscope objectives that offer exceptional image flatness and chromatic correction over 400-700nm. These objectives were designed for a 1.1" sensor to offer a large field of view and reliable image quality for high-throughput imaging applications. Ideal for life science and machine vision applications such as microscopy, flow cytometry, pharmaceutical inspection, and assembly line or fault inspection.

These objectives won the [2023 Bronze VSD Innovators Award](#).

Technical Information

Compact Tube Lens System



Standard Tube Lens System

