

Tube Lens



Stock **#18-505** **4 In Stock**

-

1

+

₹46,800

ADD TO CART

Volume Pricing	
Qty 1-4	₹46,800 each
Qty 5+	₹44,460 each
Need More?	Request Quote

Product Downloads

General

Note:

One required between camera tube and tube lens spacer/filter cube

Physical & Mechanical Properties

Length (mm):

16.15

Maximum Diameter (mm):


Optical Properties	
120.00	Focal Length FL (mm):
1X	Magnification:
Threading & Mounting	
M29 x0.75 & Ø30mm Dovetail	Mounting Threads:
Regulatory Compliance	
View	Certificate of Conformance:
China	Country of Origin:
Edmund Optics India Private Limited	Imported By:

Product Details

- Up to 43% Smaller than Industry Standard Systems
 - Optimized Performance up to 1.1" Sensor
 - Ideal for Life Science and Industrial Applications
 - Winner of the [2023 Bronze VSD Innovators Award](#)
- Edmund Optics® 120i Plan APO Infinity Corrected Objectives are infinite conjugate microscope objectives that offer exception image flatness and chromatic correction over 400-700nm. These objectives were designed for a 1.1” sensor to offer a large field of view and reliable image quality for high-throughput imaging applications. Ideal for life science and machine vision applications such as microscopy, flow cytometry, pharmaceutical inspection, and assembly line or fault inspection.
- These objectives won the [2023 Bronze VSD Innovators Award](#).

Technical Information

Compact Tube Lens System



Standard Tube Lens System

